

Detection of blood volume changes utilizing Peripheral Venous Pressure (PVP) waveform analysis



Abubakr ElSobky, Mohamed Elgamal, Ahmad Ibrahim, Kirk Shelley, M.D., PhD, Aymen Alian, M.D.

Department of Anesthesiology, Yale University School of Medicine, New Haven, CT 06520

Introduction

- Passive leg raise (PLR) is commonly used to evaluate preload fluid responsiveness and may be considered more reproducible than the fluid bolus technique¹.
- The PLR test has the same effect as mobilization of 300 cc of blood from the lower body into the cardiac chambers.
- The peripheral venous catheter is the most commonly used method of vascular access, and the peripheral venous pressure (PVP) reflects 'downstream' pressure to the right atrium. In this experiment,
- PVP waveform analysis was utilized to assess blood volume status and fluid responsiveness during mild blood volume changes, autotransfusion of 300 cc of blood (PLR test) as shown in (figure 1-A) and the loss of 600 cc of blood at -30 mmHg of Lower Body Negative Pressure (LBNP), mild hypovolemia.




Fig 1:A) PVP waveforms during LRT B) Frequency analysis of PVP waveforms

Methodology

- EKG, BP monitor, and NICOM (Cheetah Medical, MA, USA) were used to measure Thoracic Fluid Content (TFC).
- PVP waveform was analyzed using time domain analysis (mean PVP value) and frequency analysis (spectrum,2K, Hamming, Amplitude Density (AD), 93.75% overlap) with LabChart 7.3.7 (ADInstruments) as shown in (figure 1-B).

Results

- According to the PVP waveforms analysis, from the baseline to the PLR test, the mean PVP and TFC significantly increased by 48% and 1% respectively, while the DC% PVP decreased significantly by 2%.
- However, from the PLR test to -30 mmHg, the mean PVP and TFC significantly decreased by 50% and 6% respectively, while the DC% PVP significantly increased by 109%, as shown in figure 2 and 3.




Fig 2: %Change of PVP DC%, Mean PVP and PVP Cardiac frequency




Fig 3: %Change of mean PVP and TFC

Conclusion

- Transient PLR testing resulted in increased preload, which led to an increase in total fluid content and was reflected as an increase in Mean PVP at the peripheral and a reduction in the PVP DC%.
- During mild hypovolemia, then PVP DC% increased, while mean PVP and TFC decreased. These results have demonstrated the potential of using the PVP waveform as a clinical tool for monitoring changes in blood volume.

References

1) Duus, N., et al. Journal of Critical Care, 30(1). 2014